

## **NWB Format Interface**

Written by Eyal Brand and Michael Levi

### **Required Actions:**

1. Statistics and Machine Learning Toolbox.
2. Before using the class make sure to add to path @NWBRecording.
3. Please install and follow the instructions of the MatNWB package installation from its website to get the most recent version.
  - The address is : <https://github.com/NeurodataWithoutBorders/matnwb>
  - The installation should be in the following path:  
NET\timeSeriesViewer\dataRecordingObjects\helperFunctions

### **Properties:**

- pathToAllRecordings.
- pathToProcessors.
- channelNamesArr- To each recording and processor that is in the recording.
- channelNumbersArr- To each recording and processor that is in the recording.
- recordingDuration\_msArr- Each recording and processor that is in the recording.
- recordings- The info path to all recordings available.
- metadataInRec- Number of raw data in recording.
- recordingsInFile- All of the recording that are in one file.
- defaultProcessorNum- Gui supports the first recording and processor only.
- defaultRecordingNum- Gui supports the first recording and processor only.
- nwbReadFile- The file is only read, no data is loaded.
- samplingFrequencyArr- Sampling freq to each recording and processor that is in the recording.
- triggerTimestamps- To each recording and processor that is in the recording, just the data info , not loaded.
- allTriggerTimestampsLoaded- To each recording and processor that is in the recording.
- fullFilename .
- timestamps-To each recording and processor that is in the recording, just the data info , not loaded
- allTimestampsLoadedto- Each recording and processor that is in the recording.
- processorsNames.
- recordingNames.
- Metadata- The actual data.
- sample\_ms- To each recording and processor that is in the recording.

## **Methods:**

function obj = NWBRecording(recordingFile)-

This function is the class constructor. It reads the NWB file and calls getFilesAndExtract to start analyzing. It also makes sure it has all the arguments needed.

function obj =getFilesAndExtract(obj)-

This function takes the file path and read the information. Afterwards it send it to extraction function to analyze.

function [V uV ,t\_ms]=getData(obj, channels, startTime\_ms, window\_ms)-

GUI function takes the first recording and the first processor in that recoding and displays its data.

function [T\_ms]=getTrigger(obj, startTime\_ms,window\_ms)-

GUI interface applies only for one recording. The time is relative to the start time of the recording start time.

function [startIndex,endIndex]=

GetElementIndex(obj,startTime\_ms,endTime\_ms,timestampsLoaded)-

timestampsLoaded can be all timestamps / trigger timestamps etc. gives the relevant indices of the array.

function obj= extractMetaData(obj, fullFileName)-

Extracting the data from the file. goes over each recording and each processor and assign the values to the obj.

function obj=

FillInCurrentRecordingData(obj,recording,recordingName,processorNum,fullFileName)-

This function is called by extractMetaData in order to fill the values for all variables needed.

function obj= GetDefaultDataForGui(obj)-

This function is called by extractMetaData in order to fill the values for all variables needed specifically for the GUI (only the first recording and the first processor).